


I-Scan 9.0 Software Release

I-SCAN 9.0: A MODERN SOFTWARE PORTFOLIO


1. Streamlined sensor configuration and set-up
2. Simplified user interface to compare and analyze data
 - Make analysis more intuitive
 - Easy to compare data & see critical parameters
 - Manage workspace with data window display options
3. Improved Analysis options
 - Easy to compare multiple recordings
 - Rapid graph creation
 - Integrated property table
4. Web-based software updates
 - Automatic notification of update download availability



Improve productivity for data collection and analysis


STREAMLINED SENSOR CONFIGURATION & SETUP

Recommended process flow for sensor step




STREAMLINED SENSOR CONFIGURATION & SETUP

- Improved efficiency and process flow for sensor set-up
 - Save time in future sessions with remembered settings
 - Define test parameters up front
 - Ensure all steps are taken to collect accurate data




NEW DISPLAY OPTIONS FOR AN ORGANIZED WORKSPACE


- Display data, graphs or chart windows in a variety of ways
 - Docked windows to maximize real estate for displaying data
 - Split screen to view many data sets at once
 - Tab windows to prioritize important data sets
 - Auto docking of windows will simplify demo process, allowing prospects to instantly see the valuable information provided by the software



IMPROVED DISPLAY CONTROLS




Key data is displayed in the details panel



Zoom in to see details in a variety of ways:


- Slide bar
- Click to zoom region
- Scroll wheel on mouse
- Move zoom region from zoom tab

ENHANCED REGION ANALYSIS FOR DETAILED INSIGHTS



- Display data, graphs or chart windows in a variety of ways
 - Quickly create analysis graphs (pop up icon on pressure window)
 - Easily modify for more complex analysis
 - Compare key values in all regions of multiple recordings in property table


ANALYSIS OF MULTIPLE PARAMETERS



- View multiple graphs simultaneously
- Improved graph creation and modification for more complex analysis

EFFICIENT SOFTWARE UPDATES

- Simple click of a button (when connected to internet) triggers download of global bug fix updates
- Users automatically informed of version update downloads
 - No need for download links or CDs to be shipped
- Future releases will allow software changes (new maps or additional features to be pushed)




How I-SCAN 9.0 IMPROVES USABILITY

- **Complex data analysis for R&D**
 - Easier to organize workspace when comparing multiple data sets
 - Quickly configure graphs for in depth analysis
 - Conveniently display data in multi-handle applications
- **Quality Control applications**
 - Set layout for matching data to “Pass” parameters
 - Details panel and property table to quickly see key metrics


Enhance data collection with more efficiency from new user interface

ORGANIZE COMPLEX DATA ANALYSIS



Windows can be pulled from software window and displayed on second screen

SIMPLIFY QUALITY CONTROL APPLICATIONS



Docked windows maximize real estate of data images

AVAILABLE TOOLS: 7 STEP-BY-STEP TUTORIAL VIDEOS

- Located on a “Training Video” section of the Tekscan Website – **all 1-minute or less**
 - How to use the details panel and pressure analysis tools
 - How to set acquisition parameters and record data
 - How to group, tab, and reorganize windows
 - How to create graphs and tables to analyze specific regions



<https://www.tekscan.com/i-scan-90-tutorials>



Tekscan™